

Abecca recommends weekly cleaning of furniture and mattresses, as keeping Vinyl, PVC & LimberShield[™] products well maintained is essential to keep durability and visual appeal. Make sure you attend to spills immediately, wipe spill using a clean, dampened cloth and if needed a mild solution of liquid detergent and warm water or an approved cleaning solution, repeat if necessary. Please ensure all wet cleaning is followed by drying using a clean dry cloth.

The following products are Safe and effective for use on Abecca®'s products:

- Whiteley V-Wipes Hospital Grade Disinfectant
- Clinell Universal Sanitising Wipes
- **Reynard Neutral Detergent Wipes**

Please follow these recommendations to get the best use out of your furniture

If you would like to use a different cleaning product, please contact Abecca® at customerservice@abeccadesign.com before use to confirm compatibility with our Vinyl, PVC & LimberShield™ thermo-fused technology, failure to do so may void your warranty.
























For warranty purposes we highly recommend keeping a weekly cleaning log as this will be asked for along with the products used to clean the furniture when any claim is made.